


QUICK REFERENCE GUIDE



This is a quick reference guide for the **Dyna-Form Air Sure System** - (refer to user manual for additional information).

Product Code MAT/SURE/200/85/25




**Power Switch
Alarm Reset**

The Power switch simply switches the mains power to the pump on and off. When the pump detects an alarm condition, this can be cancelled by switching the pump off and then back on again.



**CPR
Valve**

Please ensure that the CPR valve is always set to normal, prior to inflating the mattress. The mattress will NOT inflate should this function not be selected. The CPR valve is only to be used in the event of a clinical emergency for priority use. However the function can also be used to deflate air rapidly from the mattress in readiness for transport. NB: The "DEFLATE" function on the control panel is also designed for this use.




LED

When the system is initially switched on or indeed should the system function be changed in any way, then the LED display will identify to the carer / user "PLEASE WAIT...SETTING PRESSURE". Please be patient and wait until the system automatically switches into the default pressure setting. This may take up to 20 minutes from the initial set up.



Ready

A green light illuminates when the system is ready for patient use. Do not use the mattress until this green light is visual.



Settings

The mode selection button is used to change the mattress cycle from "STATIC" to "ALTERNATING". The "STATIC" mode must be selected prior to transferring the patient or for when a clinical procedure requires the cells to remain stable. "ALTERNATING" mode is to be selected at all other times. A safety over-ride function automatically switches from "STATIC" to "DYNAMIC" mode after a period of 20 minutes, just in case the carer has forgotten to switch back to the required function. "COMFORT" mode allows a selection of "SOFT, MEDIUM or FIRM" accordingly. This function is simply available for patient or clinical preference. It is advised to have the unit on "Firm or Medium" mode for the heavier client (222kg / 35 stone+) "SOFT" mode can be used for hypersensitive patients or other clinical conditions deemed suitable. The "MAXIMUM INFLATION" mode can be used for a rapid increase of cell pressure for example from "SOFT to FIRM" mode.



Alarm

Other than an audible alarm, for those carers / patients that may be hard of hearing, a flashing green light also indicates that the system requires attention. During an alarm condition, then this can be silenced by pressing the reset alarm button. Once reset by the power switch above, then the alarm light will extinguish and the audible alarm will silence.



Controls

The "RESET ALARM" button is self evident for usage and simply resets the audible alarm. It is advised after using this function to switch the power off and back on again. This will reset the system and clear any non serious problem. In order to deflate the mattress completely, other than via the CPR valve, then hold down the "DEFLATE & RESET ALARM" buttons for a period of 10 seconds. This will then start deflating the mattress fully.

NB: Please ensure that all securing straps on the base of the mattress are secured onto the NON MOVING PARTS of the bed frame.

T R O U B L E S H O O T I N G



This is a troubleshooting guide for the **Dyna-Form Air Sure System** in the event of a malfunction.(refer to user manual for additional information).

Symptoms	Problems / Cause	Points to check
Constant "PLEASE WAIT... SETTING PRESSURE" LED within window.	Mattress not inflating.	Check CPR valve to normal.
"MAINS FAILURE" LED within window.	Mattress not functioning properly / deflating.	Ensure power switch is on at the socket. Please check that the power cable is not crimped / cut. If so call the helpline.
Low Pressure.	The mattress is set to mode SOFT. The CPR has been knocked / moved. There may be a leak in the system.	Change the Mode button to Medium or Firm as required. Switch CPR valve to NORMAL. Check all tubing is not kinked within the mattress. Ensure that the four connectors on the pump unit are fully locked home to the mattress.
High Pressure.	The mattress is excessively firm on a constant basis.	Set mattress to MEDIUM or SOFT mode. Check all tubing is not kinked within the mattress.
Mattress feels as though they are deflated under the patient.	Mattress is operating a natural cycle of "1 in 3". The system naturally deflates one cell in three and may be mistaken for a fault. Mattress has one or more cell failures.	Please press the static mode if unsure. This will bring all cells uniformly into a solid state, determining that the cells are functioning correctly. Switch back to required comfort setting on completion. Once ensuring that a cell failure is evident (please conduct above action first). Then call the Helpline if this persists.



Telephone Helpline / Service and Maintenance 0845 459 9831
Please document product code - MAT/SURE/200/85/25 and serial number on rear of the pump prior to calling.

Direct Healthcare Services Ltd Units 8/9, Withey Court, Western Industrial Estate, Caerphilly, South Wales, CF83 1BF, UK.

Sales Office Tel: 0845 459 9831 Fax: 0845 459 9832 E: info@direct-healthcare.org.uk W: www.direct-healthcare.org.uk

" D e l i v e r i n g t h e P r o m i s e "